

Made in
Germany

Joline®

Stop'n GO

the newest kyphoplasty
balloon catheter



twin technology®

the Revolution in Kyphoplasty



Both balloons can be individually inflated. With the individual pressure and volume adjustment, the vertebral fracture can be separately elevated - step by step. This ensures more control and safety.


Stop'n GO

Fractured structures within the vertebrae can be restored specifically, due to the individual volume and pressure control of both balloons. Step by Step.


Compression fracture with end plate transversal fracture 6 months



Controlled front side fracture of the cover plate elevated




The middle of the fracture is precisely adjusted to the geometry of the vertebrae



With both balloons, **step by step** the geometry is regained



Proximal balloon stabilizes during cementing of vertebrae



Vertebrae geometry restored and stabilized

For further information please contact us.

Joline GmbH & Co. KG
Neue Rottenburger Strasse 50
D-72379 Hechingen
Germany

T. + 49 (0) 7471 9881 0
F. + 49 (0) 7471 9881 222
E. info@joline.de
www.joline.de