







# BoneScalpel Tip Selection Guide



Tip	Part #	Approach	Purpose	Procedures
 10mm Blade, Blunt .5 x 5 mm	<b>MXB-10</b> Working Length* = 34 mm	Open	<b>Osteotomy</b> <b>En-Bloc Resection</b> <b>Bone Cutting</b>	<b>Laminoplasty</b> <b>Laminectomy</b> <b>Facetectomy</b> <b>PSO</b> <b>TLIF/PLIF</b> <b>Corpectomy</b> <b>Laminotomy</b> <b>Foraminotomy</b> <b>ACDF</b>
	<b>MXB-10LS</b> Working Length* = 145 mm	MIS		
	<b>MXB-10LC</b> Working Length* = 145 mm	MIS		
 20mm Blade, Blunt 1.0 x 6.4 mm	<b>MXB-20</b> Working Length* = 40 mm	Open		
	<b>MXB-20LC</b> Working Length* = 145 mm	MIS		
	<b>MXB-MIS-20LCRS</b> Working Length* = 100 mm	Mini-Open/MIS		
 25mm Blade, Blunt 1.0 x 6.4 mm	<b>MXB-25</b> Working Length* = 45 mm	Open		
 20mm Blade, Serrated 1.0 x 6.4 mm	<b>MXB-B1</b> Working Length* = 40 mm	Open		
 Micro Hook Shaver 1.3 x 1.8 mm	<b>MXB-S1</b> Working Length* = 57 mm	Open		
 Diamond Shaver 4.4 mm	<b>MXB-S3</b> Working Length* = 34 mm	Open	<b>360° Bone Shaving</b> <b>&amp; Bone Shaping</b> <b>Decompression</b> <b>Gross Resection</b>	<b>Corpectomy</b> <b>VCR</b> <b>Decortication</b> <b>Laminotomy</b>

\*Approximate distance from typical index finger placement to end of probe using BoneScalpel Handpeice(BCM-HP)

*If you would like further information or would like to evaluate BoneScalpel, please contact us at +1.631.694.9555*



1938 NEW HWY, FARMINGDALE, N.Y. 11735 | +1.631.694.9555 +1.631.694.3285 FAX  
 MISONIX, INC. | NASDAQ SYMBOL. MSON | MISONIX.COM

